

G136 Cavitron Plus Scaler with Tap-On Technology
Service Parts List

QUANTITY	PART NUMBER	DESCRIPTION
1	81651	Power Drive PCB Bracket Assembly
1	81635	Protective Cover, Power Entry
1	81665	Baffle (Electrical Protection)
1	81661-1	24V Power Supply Assembly
1	81646	Water Chassis Assembly
1	81843	Decorated Base
1	81632	Handpiece Holder
1	565179001	Rocker Switch, Snap-In
4	592051001	Rubber Feet, Self Adhesive
2	592055001	D-Section Seal, 1 Inch
1	568270002	Power Entry Module
1	608020001	Air Cap, Vinyl
4	586172005	#6-19 x 1/2" Hi-Lo Torx, Pan Head Screw
1	745-06002	#6 Lockwasher (For Ground Terminal)
1	586172003	#8-18 x 1/2" Hi-Lo Torx, Pan Head Screw
1	586174003	#6-19 x 3/8 Hi-Lo PPH Screw
1	62857	#6-32 x 3/8" PH Pan Head Screw
2	557034018	Fuse, 1.25A, 250V, Slo-Blo, 0.8"
1	81645	4" Ribbon Cable
1	603034001	Ferrite Block, Snap-On (For HP Cable)
1	555049007	Cable Tie (For Ferrite Block)
1	81686	DC Wire Harness
1	81662	Fuse Label
1	90176	Water Supply Label
1	629205003	Water Regulator/Solenoid (Burkert)
1	623072001	Straight Union Fitting (John Guest)
1	712-00004	Plastic Cable Tie
1	81647	Mounting Bracket
2	586174003	#6-19 x 3/8" Hi-Lo PPH Screw
1	80502-2	Knob, Molded
1	81638	Power Switch Harness
1	60366	Cable Tie (For Pot-Switch Harness)
1	81845	Cover, Decorated

G136 Cavitron Plus Scaler with Tap-On Technology
Service Parts List

6	588083001	Washers, Cup
1	81855	Molded Overlay & Label, Info Center
14	586172002	#4-24 x 1/2" Hi-Lo Screws
6	62970	Washers, 1/4" #4
1	81852	Info Center PC Board
1	592059001	Foam Seal (For "ON/OFF" Light Pipe)
2	592055001	D-Section Seal (For "ON/OFF" Light Pipe)
4	81636-2	Leveling Washers, Info Center
1	81965-2	Main Controller PC Board
1	81636-1	Leveling Washers, Main Controller
1	81645	4" Ribbon Cable
1	592056001	Shaft Seal/Nut
1	81872	Foot Pedal, Tap-On (Packed w/Synchronizing Instructions and 4 "AA" Batteries)
1	81680	Battery Door and Screw Kit
1	81681	Foot Pedal Pad (Bottom)
1	81663-2	Auxillary Foot Pedal Cable with Ferrite
6	586172005	#6-19 x 1/2" Hi-Lo Torx, Pan
1	78688	Steri-Mate Handpiece, Black (1-Pack)
1	78703	Steri-Mate Handpiece, Black (3-Pack)
1	81858	G136 Handpiece Cable, Gray
1	603034001	Ferrite Block, Snap-On (For HP Cable)
1	555049007	Cable Tie (For Ferrite Block)
1	8169201	Steri-Mate Grip 3-Pack Blue
1	8169202	Steri-Mate Grip 3-Pack Fuchsia
1	8169203	Steri-Mate Grip 3-Pack Violet
1	90125	Water Filter & Hose Assembly
1	90158	Water Filter (10 Pack)
1	81648-2	Wire Harness, Switch to Power Supply (Black Wire, White Wire)
1	81649-2	Wire Harness, Power Entry to Switch (Black Wire, White Wire, 2 Ground Wires and Two Ferrite Beads)
1	554102003	Power Cord, 115 Volt

G136 Cavitron Plus Scaler with Tap-On Technology
Service Parts List

1	78118	Power Cord, 200 to 240 Volt
1	554109001	Power Cord, 240/230 Volt-UK
1	81864	Instruction Manual (DFU) G-136
1	81865	Service Manual, G-136
1	81870	Installation Quick Guide
20120826		