

G137 Cavitron JET Plus with Tap-On Technology
Service Parts List

QUANTITY	PART NUMBER	DESCRIPTION
1	81651	Power Drive PCB Bracket Assembly
1	81635	Protective Cover, Power Entry
1	81665	Baffle (Electrical Protection)
1	81661-1	24V Power Supply Assembly
1	81713	Air / Water Chassis Assembly
1	81843	Decorated Base
1	81878	Handpiece Holder
1	565179001	Rocker Switch, Snap-In
4	592051001	Rubber Feet, Self Adhesive
2	592055001	D-Section Seal, 1 Inch
1	568270002	Power Entry Module
4	586172005	#6-19 x 1/2" Hi-Lo Torx, Pan Head Screw
1	745-06002	#6 Lockwasher (For Ground Terminal)
1	586172003	#8-18 x 1/2" Hi-Lo Torx, Pan Head Screw
1	586174003	#6-19 x 3/8 Hi-Lo PPH Screw
1	62857	#6-32 x 3/8" PH Pan Head Screw
2	557034018	Fuse, 1.25A, 250V, Slo-Blo, 0.8"
1	81645	4" Ribbon Cable
1	603034001	Ferrite Block, Snap-On (For HP Cable)
1	555049007	Cable Tie (For Ferrite Block)
1	81686	DC Wire Harness
1	81662	Fuse Label
1	81682	Air / Water Label
1	629205003	Water Regulator/Solenoid (Burkert)
1	623072001	Straight Union Fitting (John Guest)
1	712-00004	Plastic Cable Tie
1	81647	Mounting Bracket
2	586174003	#6-19 x 3/8" Hi-Lo PPH Screw
1	61679	Straight Fitting, 1/16 Barb x 10-32
1	60386	Washer, Buna & Fibre
1	81730	Pinch Plate
1	26109	Dynamic Seal
1	26105	Piston

G137 Cavitron JET Plus with Tap-On Technology
Service Parts List

1	80688	Pinch Cylinder
4	586009169	#6-32 x 7/8" Phillips Screw
1	81714	Air Manifold Assembly
1	60366	Wire Tie (For Solenoid Wires)
2	586133181	6-32 x 5/16" Phillips Screw
2	60367	Filter Assembly
1	60245	"Y" Connector, POM w/Brass Fittings
1	586161002	Pressure Relief Screw, 4-40 x 1/8"
1	81891	Powder Bowl Assembly
1	80502-2	Knob, Molded
1	81638	Power Switch Harness
1	60366	Cable Tie (For Pot-Switch Harness)
1	81879	Cover, Decorated
6	588083001	Washers, Cup
1	81881	Molded Overlay & Label, Info Center
14	586172002	#4-24 x 1/2" Hi-Lo Screws
6	62970	Washers, 1/4" #4
1	81852	Info Center PC Board
1	81720	Protective Cover (Info PC Board)
1	592059001	Foam Seal (For "ON/OFF" Light Pipe)
2	592055001	D-Section Seal (For "ON/OFF" Light Pipe)
4	81636-2	Leveling Washers, Info Center
1	81848	Main Controller PC Board
4	81636-1	Leveling Washers, Main Controller
1	81645	4" Ribbon Cable
1	592056001	Shaft Seal/Nut
1	80677	Pointer
1	628052002	O-Ring
1	80676	Powder Chamber Deflector
	80675-1	Powder Bowl Cap (White Cap only)
1	81894	Power Bowl Cap Assembly
1	81634	Bowl Mount Skirt
2	81710	L-Shape Brackets
2	60898	#6-32 x 3/8" Machine Screws

G137 Cavitron JET Plus with Tap-On Technology
Service Parts List

1	81712	Powder Bowl Housing
2	586172001	4-24 x 5/16" Hi-Lo Torx PH Screw
1	60382	Filter Element
1	60384	O-Ring
1	623071001	1/16" Elbow Barb
1	80668	Bushing
1	80672	Powder Bent Tube
1	60374	Powder Bowl
1	81711	"L" Nozzle
1	60379	Knurl Nut
1	628052001	O-Ring
1	80675	Powder Bowl Cap Assembly
1	81891	Powder Bowl Assembly, Complete
1	81872	Foot Pedal, Tap-On (Packed w/Synchronizing Instructions and 4 "AA" Batteries)
1	81680	Battery Door and Screw Kit
1	81681	Foot Pedal Pad (Bottom)
1	81663-2	Auxillary Foot Pedal Cable with Ferrite
6	586172005	#6-19 x 1/2" Hi-Lo Torx, Pan (Base to Cover)
1	8170201	JET-Mate Handpiece (1-Pack)
1	8170203	JET-Mate Handpiece (3-Pack)
1	81717	Nozzle Grip, JET-Mate Handpiece
1	81884	G137 Handpiece Cable, Black
1	603034001	Ferrite Block, Snap-On (For HP Cable)
1	555049007	Cable Tie (For Ferrite Block)
1	81075	Cleaning Wire, Jet-Mate
1	63740	Jet Air Polishing Nozzle
1	776030199	Cleaning Wire, Nozzle
1	90088	Air Filter with Fittings
1	81721	Air Filter & Hose Assembly
1	78119	Powder Jar
1	90125	Water Filter & Hose Assembly
1	90158	Water Filter (10 Pack)
2'	61693	Clear Tubing, per Foot

G137 Cavitron JET Plus with Tap-On Technology
Service Parts List

2'	61631	Red Tubing, per Foot
1'	625036010	Gray Tubing, per Foot
2'	625036012	Yellow Tubing, per Foot
2	60367	Duckbill Filter Assemblies
1	60245	"Y" Connector, POM w/Brass Fittings
4	60366	Nylon Cable Ties
1	81648-2	Wire Harness, Switch to Power Supply (Black Wire, White Wire)
1	81649-2	Wire Harness, Power Entry to Switch (Black Wire, White Wire, 2 Ground Wires and Two Ferrite Beads)
1	554102003	Power Cord, 115 Volt
1	78118	Power Cord, 200 to 240 Volt
1	554109001	Power Cord, 240/230 Volt-UK
1	81888	Instruction Manual (DFU) G-137
1	81889	Service Manual, G-137
1	81870	Installation Quick Guide
20120826		